

## Product datasheet for **MR203642L3V**

### Tgif1 (NM\_009372) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tgif1 (NM_009372) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tgif1
Synonyms:	AA959811; AI462167; Tgif
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009372
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203642).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009372.3</a> , <a href="#">NP_033398.2</a>
RefSeq Size:	1634 bp
RefSeq ORF:	819 bp
Locus ID:	21815
UniProt ID:	<a href="#">P70284</a>
Cytogenetics:	17 E1.3



[View online »](#)

**Gene Summary:**

Binds to a retinoid X receptor (RXR) responsive element from the cellular retinol-binding protein II promoter (CRBP-II-RXRE). Inhibits the 9-cis-retinoic acid-dependent RXR alpha transcription activation of the retinoic acid responsive element. May participate in the transmission of nuclear signals during development and in the adult, as illustrated by the down-modulation of the RXR alpha activities (By similarity).[UniProtKB/Swiss-Prot Function]