

Product datasheet for **MR203647A1V**

Mouse Dnajc27 (NM_153082) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Dnajc27 (NM_153082) AAV Particle
Tag:	Myc-DDK
Symbol:	Dnajc27
Synonyms:	AI639580; C330021A05Rik; Rabj; Rbj
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203647 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGACCAACGTGCCGAAGAGGAAGGAGCCCGCCAAGTCCCTCCGCATCAAAGTCATCTCCATGGGCA
ACGCCGAAGTGGGGAAAAGCTGTATTATAAAGCGCTACTGTGAGAAGAGATTTGTGTCTAAATACCTGGC
AACCATCGGGATTGACTATGGAGTCACAAAGGTCCAAGTCAGAGACAGAGAGATCAAAGTTAACATCTTT
GACATGGCTGGACATCCCTTCTCTTTGAGGTCCGAAATGAGTTTACAAGGACACACAAGGTGTGATAC
TAGTCTATGATGTTGGACAGAAAGACTCATTTGATGCCCTTGATTCTGGCTGGCAGAGATGAAGCAGGA
GCTTGGGCCTCACGGGAACATGGACAACATTGTCTTTGTTGTTTGTGCCAACAAAGATCGACTGCTCCAAA
CACCGCTGCATCGATGAAAGTGAGGGCGTCTCTGGGCTGAAAGCAAAGGTTCTGTACTTTGAAACAT
CAGCACAACTGGAGAGGGAATCAACGAGATGTTCCAGACCTTTTACTTGTCCATAGTTGACTTGTGTGA
AAATGGCGGGAAGCGTCTACTGCAAGTAGCAGTGCTAGTTTCACCAAAGAACAAGCAGACACTATTCGC
AGAATCCGAAATAGTAAGGACAGCTGGGAGATGCTAGGAGTCAGACCTGGAGCCTCAAGGGAGGAAGTTA
ATAAGGCCTACCGGAAGCTGGCTGTGCTTACCCGGACAAGTGTGTAGCACCAGGCAGTGAAGATGC
CTTCAAAGCAGTTGTGAATGCCCGACAGCCCTCTGAAAAACATCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203647 protein sequence Red=Cloning site Green=Tags(s) METNVPKRKEPAKSLRIKVISMGNAEVGKSCIIRYCEKRFVSKYLATIGIDYGVTKVQVRDREIKVNIF DMAGHPFFFEVRNEFYKDTQGVILVYDVGQKDSFDALDSWLAEMKQELGPHGNMDNIVFVVCANKIDCSK HRCIDSEGRWAESEKGFYFETSAQTGEGINEMFQTFYLSIVDLCENGGRPTASSASFTKEQADTIR RIRNSKDSWEMLGVRPGASREEVNKAYRKLAVLLHPDKCVAPGSEDAFKAVVNARTALLKNIK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_153082
ORF Size:	819 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_153082.2, NP_694722.2</u>
RefSeq Size:	4713 bp
RefSeq ORF:	822 bp
Locus ID:	217378
UniProt ID:	<u>Q8CFP6</u>
Cytogenetics:	12 A1.1
MW:	30.8 kDa