

Product datasheet for **MR203649L3V**

6720463M24Rik (BC052183) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	6720463M24Rik (BC052183) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	6720463M24Rik
Synonyms:	6720463M24Rik; AI317232
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC052183
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203649).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC052183 , AAH52183
RefSeq Size:	3866 bp
RefSeq ORF:	821 bp
Locus ID:	77744
Cytogenetics:	14 E2.2
Gene Summary:	Required for the activation of AURKA at the onset of mitosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)