

## Product datasheet for **MR203661L3V**

### **Psmb10 (NM\_013640) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Psmb10 (NM_013640) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Psmb10
Synonyms:	Mecl-1; Mecl1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013640
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203661).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013640.1</a> , <a href="#">NP_038668.1</a> , <a href="#">NP_038668.2</a>
RefSeq Size:	1270 bp
RefSeq ORF:	822 bp
Locus ID:	19171
UniProt ID:	<a href="#">O35955</a>
Cytogenetics:	8 53.06 cM



[View online »](#)

**Gene Summary:**

The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity. This subunit is involved in antigen processing to generate class I binding peptides. Plays a role in determining the T-cell repertoire for an antiviral T-cell response.[UniProtKB/Swiss-Prot Function]