

Product datasheet for **MR203662L3V**

Tpsab1 (NM_031187) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpsab1 (NM_031187) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tpsab1
Synonyms:	Mcp-7; Mcp7; Mcpt7; MMCP-7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031187
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031187.2 , NP_112464.2
RefSeq Size:	1191 bp
RefSeq ORF:	822 bp
Locus ID:	100503895
UniProt ID:	Q02844
Cytogenetics:	17 12.53 cM
Gene Summary:	Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type. May play a role in innate immunity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)