

## Product datasheet for **MR203666L3V**

### Tef (BC017689) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tef (BC017689) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tef
Synonyms:	2310028D20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC017689
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203666).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC017689</a> , <a href="#">AAH17689</a>
RefSeq Size:	4018 bp
RefSeq ORF:	821 bp
Locus ID:	21685
Cytogenetics:	15 38.25 cM
Gene Summary:	Transcription factor that binds to and transactivates the TSHB promoter. Binds to a minimal DNA-binding sequence 5'-[TC][AG][AG]TTA[TC][AG]-3' (By similarity). Also activates the telokin promoter in smooth muscle-specific and calcium-dependent manner.[UniProtKB/Swiss-Prot Function]



[View online »](#)