

## Product datasheet for **MR203676L4V**

### 6530403A03Rik (BC021443) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	6530403A03Rik (BC021443) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	6530403A03Rik
Synonyms:	1110050F08Rik; 6530403A03Rik; C6orf151
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC021443
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203676).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC021443</a> , <a href="#">AAH21443</a>
RefSeq Size:	987 bp
RefSeq ORF:	824 bp
Locus ID:	67797
Cytogenetics:	13 A3.3
Gene Summary:	Likely involved in U12-type 5' splice site recognition.[UniProtKB/Swiss-Prot Function]



[View online »](#)