

Product datasheet for **MR203686L3V**

Nubp2 (NM_011956) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nubp2 (NM_011956) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nubp2
Synonyms:	D17Wsu11e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011956
ORF Size:	828 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203686).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011956.2
RefSeq Size:	1448 bp
RefSeq ORF:	828 bp
Locus ID:	26426
UniProt ID:	Q9R061
Cytogenetics:	17 12.53 cM



[View online »](#)

Gene Summary:

Component of the cytosolic iron-sulfur (Fe/S) protein assembly (CIA) machinery. Required for maturation of extramitochondrial Fe-S proteins. The NUBP1-NUBP2 heterotetramer forms a Fe-S scaffold complex, mediating the de novo assembly of an Fe-S cluster and its transfer to target apoproteins. Negatively regulates cilium formation and structure (PubMed:23807208). [UniProtKB/Swiss-Prot Function]