

Product datasheet for MR203710L4

Plpp2 (NM_015817) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Plpp2 (NM_015817) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Plpp2
Synonyms:	L; Lpp2; PAP2-G; PAP2-gamma; PAP2c; Pp; Ppap2c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203710).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_015817
ORF Size:	828 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_015817.2 , NP_056632.2
RefSeq Size:	1646 bp
RefSeq ORF:	831 bp
Locus ID:	50784
UniProt ID:	Q9DAX2
Cytogenetics:	10 39.72 cM
Gene Summary:	The protein encoded by this gene is a lipid phosphate phosphohydrolase. It is an integral membrane protein that catalyzes the conversion of phosphatidic acid to diacylglycerol and inorganic phosphate. The transcript is expressed at high levels in lung, liver, and kidney and at low levels in brain and heart. Null mutant mice are viable and fertile and display no overt phenotypic defects. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]

Product images:



Circular map for MR203710L4