

## Product datasheet for **MR203717L3V**

### Nit2 (NM\_023175) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nit2 (NM_023175) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nit2
Synonyms:	1190017B19Rik; D16Ertd502e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023175
ORF Size:	831 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203717).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_023175.1</a> , <a href="#">NP_075664.1</a>
RefSeq Size:	1292 bp
RefSeq ORF:	831 bp
Locus ID:	52633
UniProt ID:	<a href="#">Q9JHW2</a>
Cytogenetics:	16 34.22 cM
Gene Summary:	Has a omega-amidase activity. The role of omega-amidase is to remove potentially toxic intermediates by converting 2-oxoglutarate and 2-oxosuccinamate to biologically useful 2-oxoglutarate and oxaloacetate, respectively.[UniProtKB/Swiss-Prot Function]



[View online »](#)