

Product datasheet for MR203734

Ccdc117 (NM_134033) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc117 (NM_134033) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ccdc117
Synonyms:	1110004K02Rik; 1700026O03Rik; AU018638; AU042822; AV173073; BC018601
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203734 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCCCTCGGCCGGCCCTTCAGCGGCCTCCCTCTGAGCGGCAGCGCGGACTTCTGCAGCCGCCAC
CGGCCTTCGCGGGCCGGCCTTCCCGCCAGGAGCCGCGGCCATGACCTGGCCCCGCGCCCGGTGTCCG
CGGCGCGCCGAGCAGCCGGGCGGAGAACAGCGCGCGGACGTGTTCCATTCACTGTAGAAAGAAACAC
AAGCGGTTGGCGGAGGATGATGAGTGTCCAGTAAGAAAGAAACGACTGACTGAAGCTGAAGCTCGGTGCAG
TTACTGATGAATGGGCTCTCGGCGCACATCAGGGAAGAGAAGGTCATGGGGTAAACACGTGCCTAGCAG
CCTTTCTATGCCAGCATGCTAGATGTGGTCTGTGAGGAAATGGATCAGACAACCTGGAGAACCACAGTGT
GAAGTCGCCCCGAAGGCGGCTGCAAGAGATTGAGGACAGGATAATTGACGAAGATGAAGAGTTGAAAGTG
ACAGAAATGTTAGCCACCTCCCAAGTCTTGTCTTTCTGACACCATGAAAACGGGTTTGAAGAGGGAATT
TGATGAAGTCTTCACAAAGAGGATGATTGAGTCTATGAGCCGTCCTTCTATGGAGCTTGTGCTCTGGAAA
CCTCTCCTGAACTCCTTCTGAGAAGCCAAAGCCATCATCTAGCCCTAAGAACTACAGGAGAGAGAGCC
AAGCAAAGCACGAGCTCCGGGCACTGCGTTCACAGAGGACTGAAGGGTTGCTGGAACCTCAGTGCAG
AGACGCGCCGCTTACCGTAGTCTGGAGGCAGCTACTAGCACAGAGGAGGAGATGGAGCTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203734 protein sequence
Red=Cloning site Green=Tags(s)

MAALGRPFSGPLPLSGSADFLQPPPAFAGRAFPFGAAGHDLAPRPGVVRGAPSSPGGRTARGRVSIIHCRKKH
 KRLAEDDECPVRKKRL TEAELGAVTDEWALGAHQGREGHGVNTCPSSL SMP SMLDVVCEEMDQTTGEPQC
 EVARRRLQEIEDRIIDEDEEVESDRNVSHLPSLVLSDTMKTGLKREFDEVFTKRMIESMSRPSMELVLWK
 PLPELLPEKPKPSSSPKNYRRESQAKHAAPGTAFQRTTEGLLEPQCADAPLYRSLEAATSTEEEMEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_134033

ORF Size: 834 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_134033.2](#)

RefSeq Size: 3418 bp

RefSeq ORF: 834 bp


Locus ID: 104479

UniProt ID: [Q6PB51](#)

Cytogenetics: 11 A1

MW: 30.4 kDa

Product images:



Circular map for MR203734