

## Product datasheet for MR203735

### Impa1 (NM\_018864) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Impa1 (NM_018864) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Impa1
Synonyms:	2610002K09Rik; 2900059K10Rik; AI325909
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203735 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGACCCTTGGCAGGAGTGCATGGACTATGCAGTAATCCTCGCGAGGCAAGCTGGAGAGATGATTC  
GTGAAGCTTTAAAAAATGAGATGGATGTCATGATTTAAAGTTCTCCAGCCGACTTGGTAACAGTTACTGA  
CCAAAAAGTTGAAAAATGCTCATGTCTTCTATAAAGGAAAAGTATCCATGTCACAGCTTCATTGGTGAA  
GAGTCTGTGGCAGCTGGGAGAAGACGGTCTTACAGAGCAGCCACGTGGGTTCATTGACCCATTGATG  
GAACGACTAACTTCGTGCATCGGTTCCCTTTGTAGCTGTTTCAATTGGCTTCCTTGTGAATAAAGAGAT  
GGAGTTTGAATTGTGTACAGCTGTGTGAAGATAAGATGTACACCGGCAGGAAAGGAAAGGTGCCTTT  
TGTAACGGTCAGAAGCTTCAGGTGTCCCAGCAGGAAGACATTACCAAGTCACTCTTGGTGACCGAGTTGG  
GCTCGTCCAGAAAGCCGAGACTTTACGGATCGTTCTCTCCAACATGGAAAAGCTGTGTTCCATCCCAT  
CCATGGAATCCGGAGTGTGGAACAGCTGCTGTTAATATGTGCCTTGTGGCAACGGGAGGAGCAGATGCC  
TATTATGAGATGGGAATCCACTGCTGGGACATGGCGGAGCTGGCATCATTGTACCCGAGGCAGGCGGAG  
TGCTCATGGATGTCACGGGTGGACCGTTCGATCTGATGTCTCGGAGAATAATTGCCGAAATAGTATAAC  
ATTAGCCAAAAGAATAGCCAAAGAAATTGAGATAATACCTTTGCAAAGAGACGACGAAAGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR203735 protein sequence  
 Red=Cloning site Green=Tags(s)

MADPWQECMDYAVILARQAGEMIREALKNEMDVMIKSSPADLVTVDQKVEKMLMSSIKEKYPCHSFIGE  
 ESVAAGEKTVFTEQPTWVIDPIDGTTNFVHRFPFVAVSIGFLVNKEMEFGIYVSCVEDKMYTGRKGGKAF  
 CNGQKLQVSQQEDITKSLLVTELGSSRPETLRIVLSNMEKLCSEPIHGIRSVGTAANMCLVATGGADA  
 YYEMGIHCWDMAGAGIIVTEAGGVLMDVTGGPFDLMSRRIIAANSITLAKRIAKEIEIIPLQRDDES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_018864

**ORF Size:** 834 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_018864.6](#)

**RefSeq Size:** 2601 bp


**RefSeq ORF:** 834 bp

**Locus ID:** 55980

**Cytogenetics:** 3 A1

**MW:** 30.4 kDa

**Product images:**



Circular map for MR203735