

Product datasheet for **MR203744L3V**

Tmx1 (NM_028339) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmx1 (NM_028339) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmx1
Synonyms:	2810425A04Rik; AI465319; Txndc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028339
ORF Size:	837 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203744).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028339.1 , NP_082615.1
RefSeq Size:	2418 bp
RefSeq ORF:	837 bp
Locus ID:	72736
UniProt ID:	Q8VBT0
Cytogenetics:	12 C2
Gene Summary:	May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions.[UniProtKB/Swiss-Prot Function]



[View online »](#)