

Product datasheet for **MR203750L3V**

Entpd6 (BC038126) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Entpd6 (BC038126) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Entpd6
Synonyms:	2700026H11Rik; Cd39l2; NTPDase-6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC038126
ORF Size:	834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203750).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC038126 , AAH38126
RefSeq Size:	1200 bp
RefSeq ORF:	836 bp
Locus ID:	12497
Cytogenetics:	2 G3



[View online »](#)

Gene Summary:

Catalyzes the hydrolysis of nucleoside triphosphates and diphosphates in a calcium- or magnesium-dependent manner. Has a strong preference for nucleoside diphosphates, preferentially hydrolyzes GDP, IDP, and UDP, with slower hydrolysis of CDP, ITP, GTP, CTP, ADP, and UTP and virtually no hydrolysis of ATP. The membrane bound form might support glycosylation reactions in the Golgi apparatus and, when released from cells, might catalyze the hydrolysis of extracellular nucleotides.[UniProtKB/Swiss-Prot Function]