

Product datasheet for **MR203765L3V**

Fbxo27 (NM_001163702) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Fbxo27
Synonyms:	E130008B10Rik; FBG5; Gm161
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163702
ORF Size:	840 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203765).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001163702.1 , NP_001157174.1
RefSeq Size:	1966 bp
RefSeq ORF:	843 bp
Locus ID:	233040
UniProt ID:	Q6DIA9
Cytogenetics:	7 B1



Gene Summary:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Able to recognize and bind complex-type oligosaccharides.[UniProtKB/Swiss-Prot Function]