

## Product datasheet for **MR203784L4V**

### Rasd1 (NM\_009026) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rasd1 (NM_009026) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rasd1
Synonyms:	Dexr; Dexras1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009026
ORF Size:	843 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203784).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009026.1</a>
RefSeq Size:	1609 bp
RefSeq ORF:	843 bp
Locus ID:	19416
UniProt ID:	<a href="#">O35626</a>
Cytogenetics:	11 B1.3


[View online »](#)

**Gene Summary:**

This gene encodes a member of the Ras superfamily of small GTPases and is induced by dexamethasone. The encoded protein is an activator of G-protein signaling and acts as a direct nucleotide exchange factor for Gi-Go proteins. This gene may play a role in dexamethasone-induced alterations in cell morphology, growth and cell-extracellular matrix interactions. [provided by RefSeq, Nov 2015]