

Product datasheet for **MR203799L3V**

Rsrc1 (BC038303) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rsrc1 (BC038303) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rsrc1
Synonyms:	1200013F24Rik; SRrp53
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC038303
ORF Size:	843 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203799).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC038303 , AAH38303
RefSeq Size:	1182 bp
RefSeq ORF:	845 bp
Locus ID:	66880
Cytogenetics:	3 E1
Gene Summary:	Plays a role in pre-mRNA splicing. Involved in both constitutive and alternative pre-mRNA splicing. May have a role in the recognition of the 3' splice site during the second step of splicing (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)