

Product datasheet for MR203845

Haghl (NM_026897) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Haghl (NM_026897) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Haghl
Synonyms:	1500017E18Rik; 2810014I23Rik; C330022E15Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203845 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGTCAAAGTCATCCCGTACTTGAAGACAACACTACATGTACCTGATCATTGAGGAACACACGCGGG
AGGCCGTAGCTATAGACGTGGCTGTGGCCGAGAGGTTGCTGGAGATTGCGGGTCGCGAGGGAGTGTCTCT
GACCATGGTGTGTCCACCCATCATCACTGGGACCACACACGTGAAATGCGGAGCTGGCACACATTCTG
CCAGGATTGGCTGTCTTGGGAGCAGATGAGCGCATTTGCGCTCTAACCCGAAGGCTGGAACATGGAGAAG
GGCTTCAGTTTGGGGCCATCCACGTGCGATGCCTCTTGACGCCAGGCCACACCTCTGGCCACATGAGCTA
CTTCTGTGGGAGGATGACTGCCCGGACTACCCAGCCCTTTCTCTGGGGATGCTCTGTCAGTGGCAGGC
TGTGGCTGGCACCTGGAGGACACTGCCAGCAGATGTACCAGAGCTTGGCCAAGACTCTGGGTACCCCTT
CACCCGAGACGAAGGTGTTCTGTGGTCATGAGCATACGCTGAGCAATCTTGAGTTTGACAGAAAAGTGGA
GCCCTGCAACGAACACGTGCAAGCTAAGCTGTCCTGGGCCAGGAGAGGGATGACGAGGATATTCCTACT
GTGCCTTCAGCGCTGGGCGAGGAGCTCATGTACAACCTTTCTGAGGGTGACAGAGGATGCAGTCCGTG
CGTTCACAGGCCAGGTGGCCCCAGCCAGGTCTGGAGGCCCTGTCAGGGAGCGAGCGCGCTTCCAGCC
TGCTGTGGAGCCACCACAGCCGAGTTTCGGGCACTCTGGCCTTACAGTGGGGGCTTCTGAGCACACAC
CAGAAAAAG

ACGCGTACGCGGCCGCTCGAGCAGAAAACACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203845 protein sequence
Red=Cloning site Green=Tags(s)

MKVKVIPVLEDNYMYLIIIEHTREAVAIDVAERLLEIAGREGVSLTMVLSTHHHWDHTRGNAELAHIL
 PGLAVLGADERICALTRRLEHGEGQLQFGAIHVRCLLTPGHTSGHMSYFLWEDDCPDSPALFSGDALSVAG
 CGWHLLEDTAQQMYQSLAKTLGTLPPETKVF CGHEHTLSNLEFAQKVEPCNEHVQAKLSWAQERDDEDIPT
 VPSTLGEELMYPFLRVTEDAVRAFTGQVAPAQVLEALCRERARFQPAVEPPQPQVRALLALQWGLLSTH
 QKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_026897

ORF Size: 852 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026897.3](#)

RefSeq Size: 1321 bp

RefSeq ORF: 852 bp

Locus ID: 68977


UniProt ID: [Q9DB32](#)

Cytogenetics: 17 A3.3

MW: 31.5 kDa

Gene Summary: Hydrolase acting on ester bonds.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203845