

Product datasheet for **MR203869**

Plet1 (NM_029639) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Plet1 (NM_029639) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Plet1
Synonyms:	0610037B23Rik; 1600029D21Rik; AA407923; AU043048
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203869 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCTGTCGCTCCGCTCCTTGCTTCCACACCTGGGACTGTTCTGTGCCTGGCTCTGCACTTATCCCCCT
 CCCTCTTGCCAGTGATAATGGTCCTGCGTGGTCCTTGATAACATCTACACCTCCGACATCTTGAAAT
 CAGCACTATGGCTAACGTCTCTGGTGGGATGTAACCTATACAGTGACGGTCCCGTGAACGATTGAGT
 AGTGCCGTGATCCTGAAAGCAGTGAAGGAGGACGACAGCCAGTGGGCACCTGGAGTGAACATATGAGA
 AGTGCAACGACAGCAGTGTCTACTATAACTTGACATCCCAAAGCCAGTCGGTCTTCCAGACAACTGGAC
 AGTTCCTACTTCCGAGGATGTGACTAAAGTCAACCTGCAGGTCCATCGTCGTAATCGCACAGCCTCA
 AAGTCATCCGTGAAAATGGAACAAGTACAACCCCTCAGCCTCAACCCCTATTCTGAGAGTTCTGAGACCA
 GCCAGACCATAAACACGACTCCAAGTGAACACAGCCAAGACTACAGCCAAGGACACAGCCAATACCA
 AGCCGTGACCACAGCCAATACCACAGCCAATACCACAGCCGTGACCACAGCCAAGACCACAGCCAAAAGC
 CTGGCCATCCGCACTCTCGGCAGCCCCCTGGCAGGTGCCCTCCATATCCTGCTTGTCTTCTCATTAGTA
 AACTCCTCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR203869 protein sequence
 Red=Cloning site Green=Tags(s)



MLSLRSLPHLGLFLCLALHLSPSLASDNGSCVVDNIYTSIDILEISTMANVSGGDVITYTVTPVNDSV
 SAVILKAVKEDDSPVGTWSGTYEKCNDSVVYNYLTSQSQSVFQTNWTPTSEDVTKVNLQVLIVVNRAS
 KSSVKMEVQPSASTPIPESETSQTINTPTVNTAKTTAKDTANTTAVTTANTTAVTTAKTTAKS
 LAIRTLGSPLAGALHILLVFLISKLLF

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 029639

ORF Size: 711 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029639.2](#), [NP_083915.2](#)

RefSeq Size: 1846 bp

RefSeq ORF: 714 bp

Locus ID: 76509


UniProt ID: [Q8VEN2](#)

Cytogenetics: 9 A5.3

MW: 25 kDa

Gene Summary: Modulates leading keratinocyte migration and cellular adhesion to matrix proteins during a wound-healing response and promotes wound repair. May play a role during trichilemmal differentiation of the hair follicle.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203869