

Product datasheet for **MR203869L3V**

Plet1 (NM_029639) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plet1 (NM_029639) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plet1
Synonyms:	0610037B23Rik; 1600029D21Rik; AA407923; AU043048
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029639
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203869).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029639.2 , NP_083915.2
RefSeq Size:	1846 bp
RefSeq ORF:	714 bp
Locus ID:	76509
UniProt ID:	Q8VEN2
Cytogenetics:	9 A5.3
Gene Summary:	Modulates leading keratinocyte migration and cellular adhesion to matrix proteins during a wound-healing response and promotes wound repair. May play a role during trichilemmal differentiation of the hair follicle.[UniProtKB/Swiss-Prot Function]


[View online »](#)