

## Product datasheet for **MR203884L4V**

### Tia1 (BC023813) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tia1 (BC023813) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tia1
Synonyms:	2310050N03Rik; AI256674; mTIA-1; TIA-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC023813
ORF Size:	855 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203884).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC023813</a> , <a href="#">AAH23813</a>
RefSeq Size:	1455 bp
RefSeq ORF:	857 bp
Locus ID:	21841
Cytogenetics:	6 D1
Gene Summary:	Involved in alternative pre-RNA splicing and regulation of mRNA translation by binding to AU-rich elements (AREs) located in mRNA 3' untranslated regions (3' UTRs). Possesses nucleolytic activity against cytotoxic lymphocyte target cells. May be involved in apoptosis (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)