

Product datasheet for **MR203907L4V**

Aspdh (NM_026690) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aspdh (NM_026690) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aspdh
Synonyms:	0610012D14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026690
ORF Size:	864 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026690.1 , NP_080966.1
RefSeq Size:	966 bp
RefSeq ORF:	864 bp
Locus ID:	68352
UniProt ID:	Q9DCQ2
Cytogenetics:	7 B3
Gene Summary:	Specifically catalyzes the NAD or NADP-dependent dehydrogenation of L-aspartate to iminoaspartate.[UniProtKB/Swiss-Prot Function]



[View online »](#)