

Product datasheet for **MR203910**

Sh3glb2 (BC028859) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sh3glb2 (BC028859) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sh3glb2
Synonyms:	SPAS-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203910 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACTCAACATGAAGAAGCTGGCGTCGGACGCGGGCATCTTCTTCACTCGGGCGGTGCAGTTCACAG
AGGAGAAGTTTGGCCAGGCTGAGAAGACGGAGCTTACGCCCACCTTGAAAACCTCCTAGCTCGGGCAGA
CAGCACCAAGAAGTGGACAGAGCGGATCCTGAGGCAGACCGAGGTGCTGCTGCAGCCCAACCCAGTGCT
CGAGTGGAGGAGTTCCTATATGAGAAGCTGGACAGAAAGGTGCCCTCGAGAGTACCAATGGGGAGCTGC
TGGCTCAGTACATGGCGGAGGCAGCCAGCGAGCTGGGCCCCAGCACTCCCTACGGGAAGACGCTGATCAA
GGTGTGAGAAGCTGAGAAGCGCCTCGGAGCAGCAGAGCGGGATTCATTACACTGCCTCCCTCAGCTTC
CTCACACCCCTGCGGAATTCCTAGAGGGGGACTGGAAAACGATTCGAAGGAGAGGCGGCTCCTGCAGA
ACCGGCGTCTTGACCTGGATGCCTGCAAAGCCCGCTAAAGAAGGCCAAGGCAGCCGAAGCCAAAGCCAC
GCTTTGGAACGATGAAGTCGACAAGGCTGAGCAGGAGCTTCGAGTGGCGCAGACAGATTTGACCGGCAG
GCAGAAGTGACCCGCTCCTGCTGGAGGGGATCAGCAGTACTCATGTGAACACCTTCGCTGCCTCCATG
AGTTCGTAAGTCTCAGACAACCTACTATGCACAGTGTACCGCCACATGCTGGATCTGCAGAAACAGCT
AGGCAGGTGTGGCTGGACCCAGGTCTCTGAGGCCCCCTATCTCCAGACCCAGTGAAGTGGCCCAATGCTT
CTGGCCAAGGGCCTTCTTTC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203910 protein sequence
 Red=Cloning site Green=Tags(s)

MDFNMKKLASDAGIFFTRAVQFTEEKFGQAEKTELDAHFENLLARADSTKNWTERILRQTEVLLQPNPSA
 RVEEFLYEKLDKRVPSRVNTGELLAQYMAEAASELGPSTPYGKTLIKVSEAEKRLGAAERDFIHTASLSF
 LTPLRNFLEGDWKTISKERRLLQNRRLDL DACKARLKKAKAAEAKATLWNDEVKAEQELRVAQTEFDRQ
 AEVTRLLLEGISSTHVNHLRCLHEFVKSQTTYAQCYRHMLDLQKQLGRCGWTQVSEAPYLPDPVSGPML
 LAKGLLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC028859

ORF Size: 861 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC028859](#), [AAH28859](#)

RefSeq Size: 1776 bp


RefSeq ORF: 863 bp

Locus ID: 227700

Cytogenetics: 2 B

MW: 32.8 kDa

Product images:



Circular map for MR203910