

Product datasheet for MR203930

Pbld (BC019997) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pbld (BC019997) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pbld
Synonyms:	0610038K03Rik; Mawbp; Pbld
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203930 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTCCCAATCTTCATAGCAGATGCGTTCCTGCGACAGCTTTCGGGGCAATCCTGCCGAGTCT
GCCTCCTGGAGCGCACATTGGAGGAAGACGCCATCAGCAAATTGCGAGGGAGATGAACCTCTCAGAAAC
GGCTTTCATCAGGAACTGCAGCCGACTGACAGCTTACACAAAGCTCTCGCTTGGACTGAGATGGTTT
ACGCCAGTGAGTGAAGTCCCCTTGTGCGGCCACGCCACTGGCCTCTGCGGCTGTGCTGTTTACAAAA
TACAGAACAGGAACAGCACCCCTGACCTTGTGACAATGAGCGGGGAGCTAAAGGCCAGAAGAGCAGAAGA
TGGGATAGTCTGGACTTTCCTGTCTACCAACCTTCCCGAGGACTTCCATGAAGTTGAAGACTTGATA
AAGGCAGCCATAGGTGACACCCTGGTCCAGGACATCCGGTACTCTACAGATACCAGAAACTCCTGGTCC
GGCTCAGCGATTCTATGACAGGTCTTCTAGAGAGCCTGAAGGTGAACACAGAGCCTCTGCCTGCAAT
TGAGAAGACAGGGAAGGTGAGAGGCCTCATCTCACTGTCAAAGGAGAGCCTGGGGGGCAGACGGCCCCG
TATGACTTCTACTCCAGATACTTCGCGCCGTGGTTGGTATTGCTGAAGACCCGGTGACAGGATCTGCAC
ACACCGTTCTCAGCAGTTACTGGTCCAGCAGCTCAGGAAGAAAGAGATGCGAGCCTTTCAGTGTCCCCG
CCGAGGAGGGGAGCTGGACATTTCCCTGCGACCTGATGGGCGCTTGACATAAAGGGCGGAGCTGTTATT
GTCTTGGAGGGCACGCTGACCGCC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203930 protein sequence
Red=Cloning site Green=Tags(s)

MKLPIFIADAFTATAFRGNPAAVCLLERTLEEDAHQQIAREMNLSETAFIRKLQPTDSFTQSSRFGLRWF
 TPVSEVPLCGHATLASAAVLFHKIQRNSTLTFVTMSGELKARRAEDGIVLDFPVYPTFPQDFHEVEDLI
 KAAIGDTLVQDIRYSTDRKLLVRLSDSYDRSFLESCLKVNTEPLPAIEKTGKVRGLILTVKGEPGGQTAP
 YDFYSRYFAPWVGIAEDPVTGSAHTVLSSYWSQQLRKKEMRAFCSSRRGGELDISLRPDGRVDIKGGA
 VLEGLTLTA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: BC019997

ORF Size: 864 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC019997](#), [AAH19997](#)

RefSeq Size: 1956 bp


RefSeq ORF: 866 bp

Locus ID: 68371

Cytogenetics: 10 B4

MW: 32 kDa

Product images:



Circular map for MR203930