

Product datasheet for **MR203931L4V**

Tra2b (NM_009186) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Tra2b (NM_009186) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Tra2b |
| Synonyms: | 5730405G21Rik; D16Ertd266e; Sfrs10; SIG-41; Silg41; TRA2beta |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_009186 |
| ORF Size: | 864 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR203931). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_009186.4 |
| RefSeq Size: | 2135 bp |
| RefSeq ORF: | 867 bp |
| Locus ID: | 20462 |
| UniProt ID: | P62996 |
| Cytogenetics: | 16 13.18 cM |



[View online »](#)

Gene Summary:

Sequence-specific RNA-binding protein which participates in the control of pre-mRNA splicing. Can either activate or suppress exon inclusion. Acts additively with RBMX to promote exon 7 inclusion of the survival motor neuron SMN2. Activates the splicing of MAPT/Tau exon 10. Alters pre-mRNA splicing patterns by antagonizing the effects of splicing regulators, like RBMX. Binds to the AG-rich SE2 domain in the SMN exon 7 RNA. Binds to pre-mRNA (By similarity).[UniProtKB/Swiss-Prot Function]