

Product datasheet for MR203933

Morf4l2 (NM_001168230) Mouse Tagged ORF Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


Product Type:	Expression Plasmids
Product Name:	Morf4l2 (NM_001168230) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Morf4l2
Synonyms:	2410017O14Rik; mKIAA0026; Mrgx; Sid393p
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	<pre>>MR203933 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)</pre>
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C

ATGAGTTCCAGAAAGCAGGCTTCTCAAACTCGTGGACAACAATCTGCTGAAGAAGACAACTTTAAGAAAC CAACTCGAAGCAATATGCAGAGAAGTAAGATGAGAGGATCTGCCTCGGGAAAGAAGTCAGCTGGTTCGCA GCCAAAGAATCTCGATCCAGCCCTGCCCGGAAGATGGGGAGGTCGCTCTGCTGAGAACCCCCCTTCCGGT TCTGTGCGGAAGACCAGGAAGAACAAGCAGAAGGCTCCTGGCAACGGAGACGGAGGCAGTACCAGTGAAG TCCCCCAGCCCCTCGGAAGAACAAGGCACGGGCTGACCCCACTGTGGAGAACGGAGGCAGTACCAGTGAAG TAGGATGGAGGTGAAGGTGAAGATCCCTGAAGAATTAAAACCGTGGCTGGTGGAGGACTGGGACTTCAAGAG TAGGATGGAGGTGAAGGTGCCCCGCTGTCCCCGCTAAAAAGAATGTCGATGCCATTCTTGAGGAGTATGCCA ATTGTAAGAAGTCGCAGGGAAATGTTGATAATAAGGAGTACGCAGTTAATGAAGTTGTAGGAGGATAAA AGAGTATTTCAATGTGATGCTGGCCACTCAGCTGCCGTGTACAAGTTTGAAAGGCCTCAGTATGCTGAGA GAATTGGGGCAATGTTGGCCCATGCCGCAGTCATGCGGCGCCCACCCTCCTGAGATTATTCGTGA GAATTGGGGCAATGTTGGCCTATACGCCCCTTGATGAGAAAAGCCTGGCATTATTGCTGGGCTATCTGCA TGGTTGGGCCAATGTTGGCCAAAGAATTCTGCCTCTCTGTTTACTGCCAGTGAGTTACAAAGTGGCTTCT GCTGACTATCATCGCAAAGCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAG**GTTTAA**



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Morf4l2 (NM_001168230) Mouse Tagged ORF Clone - MR203933


This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

ORIGENE Morf4l2 (NM_001168230) Mouse Tagged ORF Clone – MR203933	
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Metho	 d: 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001168230.1, NP 001161702.1</u>
RefSeq Size:	1741 bp
RefSeq ORF:	867 bp
Locus ID:	56397
UniProt ID:	<u>Q9R0Q4</u>
Cytogenetics:	X F1
MW:	32.2 kDa
Gene Summary:	Component of the NuA4 histone acetyltransferase complex which is involved in transcriptional activation of select genes principally by acetylation of nucleosomal histone H4 and H2A. This modification may both alter nucleosome - DNA interactions and promote interaction of the modified histones with other proteins which positively regulate transcription. This complex may be required for the activation of transcriptional programs associated with oncogene and proto-oncogene mediated growth induction, tumor suppressor mediated growth arrest and replicative senescence, apoptosis, and DNA repair. The NuA4 complex ATPase and helicase activities seem to be, at least in part, contributed by the

complex ATPase and helicase activities seem to be, at least in part, contributed by the association of RUVBL1 and RUVBL2 with EP400. NuA4 may also play a direct role in DNA repair when directly recruited to sites of DNA damage. Also component of the MSIN3A complex which acts to repress transcription by deacetylation of nucleosomal histones (By similarity).[UniProtKB/Swiss-Prot Function]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Circular map for MR203933

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US