

Product datasheet for **MR203946L3V**

Yipf4 (BC003936) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Yipf4 (BC003936) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Yipf4
Synonyms:	2310034L04Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC003936
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203946).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC003936 , AAH03936
RefSeq Size:	2007 bp
RefSeq ORF:	869 bp
Locus ID:	67864
Cytogenetics:	17 E2
Gene Summary:	Involved in the maintenance of the Golgi structure.[UniProtKB/Swiss-Prot Function]



[View online »](#)