

Product datasheet for **MR203948L4V**

Pnp (NM_013632) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pnp (NM_013632) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pnp
Synonyms:	AL024301; AU015798; Np; Np-1; Np-2; Pnp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013632
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013632.4 , NP_038660.1
RefSeq Size:	2941 bp
RefSeq ORF:	870 bp
Locus ID:	18950
UniProt ID:	P23492
Cytogenetics:	14 26.31 cM
Gene Summary:	The purine nucleoside phosphorylases catalyze the phosphorolytic breakdown of the N-glycosidic bond in the beta-(deoxy)ribonucleoside molecules, with the formation of the corresponding free purine bases and pentose-1-phosphate.[UniProtKB/Swiss-Prot Function]



[View online »](#)