

Product datasheet for **MR203954L3V**

Pih1d1 (NM_029406) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pih1d1 (NM_029406) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pih1d1
Synonyms:	1110061L23Rik; 4933413A04Rik; Nop17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029406
ORF Size:	873 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203954).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029406.1 , NP_083682.1
RefSeq Size:	1196 bp
RefSeq ORF:	873 bp
Locus ID:	68845
UniProt ID:	Q9CQJ2
Cytogenetics:	7 B3



[View online »](#)

Gene Summary:

Involved in the assembly of C/D box small nucleolar ribonucleoprotein (snoRNP) particles. Recruits the SWI/SNF complex to the core promoter of rRNA genes and enhances pre-rRNA transcription. Mediates interaction of TELO2 with the R2TP complex which is necessary for the stability of MTOR and SMG1. Positively regulates the assembly and activity of the mTORC1 complex.[UniProtKB/Swiss-Prot Function]