

Product datasheet for **MR203974L3V**

Mff (NM_029409) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mff (NM_029409) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mff
Synonyms:	5230400G24Rik; AI314724
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029409
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203974).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029409.2 , NP_083685.2
RefSeq Size:	1920 bp
RefSeq ORF:	876 bp
Locus ID:	75734
UniProt ID:	Q6PCP5
Cytogenetics:	1 C5



[View online »](#)

Gene Summary:

Plays a role in mitochondrial and peroxisomal fission. Promotes the recruitment and association of the fission mediator dynamin-related protein 1 (DNM1L) to the mitochondrial surface. May be involved in regulation of synaptic vesicle membrane dynamics by recruitment of DNM1L to clathrin-containing vesicles.[UniProtKB/Swiss-Prot Function]