

## Product datasheet for MR203985L4

### Tmem41b (NM\_153525) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem41b (NM_153525) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Tmem41b
Synonyms:	1500015G02Rik; 1500031M19Rik; AU018901; D7Erttd70e; D7Erttd743e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203985).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_153525
ORF Size:	876 bp



[View online »](#)

