

## Product datasheet for **MR203991L4V**

### 2310037124Rik (BC092275) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                                  |                                   |
|----------------------------------|-----------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles              |
| <b>Symbol:</b>                   | 2310037124Rik                     |
| <b>Synonyms:</b>                 | 2310037124Rik; B930009D17         |
| <b>Mammalian Cell Selection:</b> | Puromycin                         |
| <b>Vector:</b>                   | pLenti-C-mGFP-P2A-Puro (PS100093) |
| <b>Tag:</b>                      | mGFP                              |
| <b>ACCN:</b>                     | BC092275                          |
| <b>ORF Size:</b>                 | 873 bp                            |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(MR203991).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**RefSeq:** [BC092275](#), [AAH92275](#)

**RefSeq Size:** 1784 bp

**RefSeq ORF:** 875 bp

**Locus ID:** 69612

**Cytogenetics:** 15 F1



**Gene Summary:**

As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.  
[UniProtKB/Swiss-Prot Function]