

Product datasheet for **MR203994L3V**

Imp4 (NM_178601) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Imp4 (NM_178601) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Imp4
Synonyms:	AA409888; AV031295; D1Wsu40e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178601
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203994).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178601.2
RefSeq Size:	2249 bp
RefSeq ORF:	876 bp
Locus ID:	27993
UniProt ID:	Q8VHZ7
Cytogenetics:	1 13.22 cM
Gene Summary:	Component of the 60-80S U3 small nucleolar ribonucleoprotein (U3 snoRNP). Required for the early cleavages during pre-18S ribosomal RNA processing (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)