

Product datasheet for MR203996

Fam84a (NM_029007) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam84a (NM_029007) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fam84a
Synonyms:	2310003N02Rik; 4731402F03Rik; AW125753; Nse1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203996 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCAACCAACTGGACCGCATCACCCACCTCAACTACAGCGAGTTGCCACTGGGGACCCATCTGGGA
TTGAGAAGGACGAGCTGCGAGTCGGGGTTGCCTACTTCTTCTCGGATGAGGAGGAGACCTGGACGAACG
CGGGCAGCCGGACAAGTTTGGTGTGAAGGGCCCCCAGGTTGCAGCCCTTGCCAGAGAGCCCCAGCCGC
CACCACCACCTGCTGCACCAGCTCGTCTCAACGAGACTCAGTTCTCCGCCTCCGGGGCCAGGAAT
GCATCTTTTCAAAGTGACCGGGGCCCTCAGGGCCGACTTGAGTGTCTACGGGTCACTGCGCTGCC
GGCGATCTGCGAGCCAGCGACCTGCTAGAGCTCTTGTGGCTACAGCCGGCGACCGAGCAGCCCGGCCA
GCCCTCACTGGGCCGTCTACGTGGGCGGGCAGATCATCCACTGCACCAAGGCGAGATCCGCCAGG
ACAGCCTGTACCAGGCGGGCGGCCAACGTGGGCCGGTGGTGAATAGCTGGTACCGCTACCGCCCGCT
GGTGGCCGAGCTAGTGGTGCAGAACGCCTGTGGCCACCTGGGCCTCAAGAGCGAGGAGATCTGCTGGACG
AACTCCGAGAGCTTCGCTGCCTGGTGGCGCTTTGGAAAGCGCGAGTTCAAGGCTGGCGGGGAGGTCCCCG
CAGGAACGCAGCCCGCAACAGCAGTATTATCTCAAAGTGCATCTGGAGGAGAACAAAGTCCACACGGC
CCGTTCCACAGCCTGGAGGACCTATCCGCGAGAAGCGTGCATAGACGCCAGTGGCCCGCTGCGAGTG
CTGCAGGAGCTTGAGACTTTGTGGACGACAAGGAG

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203996 protein sequence
Red=Cloning site Green=Tags(s)

MGNQLDRITHLNYSELPDPSGIEKDELRVGVAYFFSDEEEDLDERGQDPKFGVKGPPGCSPCPESPSR
 HHHHLLHQLVLNETQFSAFRGQECIFSKVTGGPQGADLSVYAVTALPAICEPGDLELLWLQPATEQPAP
 APHWAVYVGGGQIIHLHQGEIRQDSL YQAGAANVGRVVNSWYRPLVAELVVQNACGHLGLKSEEICWT
 NSESF AAWCRFGKREFKAGGEVPAGTQPPQQYYLKVHLEENKVHTARFHSLEDLIREKRRIDASGRLRV
 LQELED FVDDKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_029007

ORF Size: 879 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029007.2](#), [NP_083283.2](#)

RefSeq Size: 3402 bp

RefSeq ORF: 879 bp

Locus ID: 105005


UniProt ID: [Q9D650](#)

Cytogenetics: 12 A1.1

MW: 32.7 kDa

Gene Summary: May play a role in cell morphology and motility.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203996