

Product datasheet for **MR203999L4V**

Sdcbp2 (NM_145535) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sdcbp2 (NM_145535) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sdcbp2
Synonyms:	MGC11704
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145535
ORF Size:	879 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145535.1 , NP_663510.1
RefSeq Size:	1412 bp
RefSeq ORF:	879 bp
Locus ID:	228765
UniProt ID:	Q99JZ0
Cytogenetics:	2 G3
Gene Summary:	Binds phosphatidylinositol 4,5-bisphosphate (PIP2). May play a role in the organization of nuclear PIP2, cell division and cell survival.[UniProtKB/Swiss-Prot Function]



[View online »](#)