

Product datasheet for MR204027

Rps2 (BC094058) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rps2 (BC094058) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rps2
Synonyms:	Llrep3; S4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204027 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGATGACGCCGGTGCAGCGGGAGGGCCCGGAGGACCCGGGGGCCAGGATTAGGTGGCCGGCGG
GCTTCCGCGGAGGATTCGGCAGCGGTCTTAGGGGCCGTGGTCTGGCCGAGGCCGTGGCCGTGGTCGAGG
CCGCGGGGCTCGTGGAGGTAAGCTGAAGACAAGGAGTGGATCCCGTCACCAAGCTGGCCGTCTGGTT
AAGGACATGAAGATCAAGTCCTTGGAGGAGATCTACCTGTTCTCCCTGCCATTAAGGAGTCTGAGATCA
TTGATTTCTTCTGGGTGCATCCCTAAAGGATGAGGTTCTGAAAATCATGCCAGTGCAGAAGCAGACTCG
GGCTGGCCAGCGGACCAGGTTCAAGGCTTTCGTCGCTATTGGGGACTACAATGGTCACGTTGGTCTTGGT
GTTAAGTGCTCCAAGGAGTTGCCACTGCCATCCGAGGGGCCATCATCTTGGCCAAGCTTCCATCGTCC
CTGTGCGGAGAGGCTACTGGGGAAACAAGATTGGCAAGCCCCACACTGTTCCATGCAAGGTGACAGGCCG
CTGTGGCTCTGTGCTGGTGCCTCATCCCTGCCCCAGAGGCACTGGCATTGTCTCTGCTCCTGTGCC
AAGAAGCTCCTGATGATGGCCGGCATAGATGACTGCTACACTTCAGCCAGAGGCTGCACTGCCACCCTGG
GCAACTTGTCTAAGGCCACCTTTGATGCCATCTCAAGACTTACAGCTACCTGACCCCGACCTCTGGAA
AGAGACTGTCTTCAACAAGTCTCCTTATCAGGAATTCACGGATCATCTTGTGAAAACCCACACCAGAGTC
TCTGTTACAGGACCCAGGCTCCAGCTGTGGCTACCACA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204027 protein sequence
Red=Cloning site Green=Tags(s)

MADDAGAAGGPGGPGGLGGRGGFRGGFGSGLRGRGRGRGRGRGRGARGGKAEDKEWIPVTKLGRLLV
 KDMKIKSLEEIYLFSLPIKESIIDFFLGASLKDEVLKIMPVQKQTRAGQRTRFKAFVAIGDYNGHVGLG
 VKCSKEVATAIRGAIILAKLSIVPVRGGYWGNIKIGKPHTVPCVKVTGRCGSVLVRLLIPAPRGTGIVSAPVP
 KKLLMMAGIDDCYTSARGCTATLGNFAKATFDIAISKYSYLTPLDWKETVFTKSPYQEFTHLVKTHTRV
 SVQRTQAPAVATT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: BC094058

ORF Size: 879 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC094058](#), [AAH94058](#)

RefSeq Size: 957 bp


RefSeq ORF: 881 bp

Locus ID: 16898

Cytogenetics: 17 A3.3

MW: 31.2 kDa

Product images:



Circular map for MR204027