

Product datasheet for **MR204044A1V**

Mouse Otud6b (BC031474) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Otud6b (BC031474) AAV Particle
Tag:	Myc-DDK
Symbol:	Otud6b
Synonyms:	2600013N14Rik; AU015433
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204044 representing BC031474 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGAGGTCGTGGCTGAAGAGCTCGATGATGAGGAGCAGCTGGTGAGAAGGCATCGGAAGGAGAAGA
AGGAGCTGCAAGCCAAAATTCAGGGAATGAAGAACGCTGTTCCCAAAAACGACAAAAGAGGAGGAAGCA
GCTCACGGAAGATGTTGCTAAATGGAAAGAGAAATGGAGCAAAAACACAGGGAAGAAGCTGGAGCAATTG
AAGCAATTGACTTTCAAGGACAGTAAGATAGATTCTGTTGCTGTTAACATTTCAAACCTTGGTACTTGAGA
ATCAACCACCTCGGATTTCAAAGCACAAAAGAGACGGGAAAAGAAGGCTGCATTGGAAAAGAGCGGGA
AGAAAGGATAGCAGAGGCTGAAATTTGAAACTTATCTGGAGCCAGACATCTTGAGAGTGA AAAACTTGCT
CAAATATTGGCAGCCAGAGAATTGAAAATTAACAGATTCCATCTGATGGTCACTGTATGTATGGAGCAC
TTGAGGATCAGCTGAGAGAACAAGACTGTGCGCTGACTGTGGCTTCCCTGAGGAGGCAGACTGCTGAGTA
CATGCAAACCCATTTCAGACGACTTCTGCGCTTCTAACAAACCCAGCACAGGAGATATGTACACGCCA
GAAGAGTTTGGGAAGTACTGTGACGATATTGTGAACACAGCGGCATGGGGCGGTCAGCTCGAGCTAAGAG
CTTTGTCTCACATCTTACAAACCAATCGAGATACTACAGGCAGATGCTCCTCCTATTATAGTCGGTGA
AGAGTATCCGAGGAATCCATTAGTACTAGTATACATGAGGCACGCGTATGGCTTAGGAGAACATTATAAT
TCTGTTACACGGTTGGTGAACCTAGCTACTGAAAATTGCAGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204044 representing BC031474 Red=Cloning site Green=Tags(s) MEEVVAEELDDEEQLVRRHRKEKELQAKIQGMKNAVPKNDKKRRKQLTEDVAKLEREMEQKHREELEQL KQLTFKDSKIDSVAVNINLVLNQPPRISKAQKRREKKAALKEKEREERIAEAEIENLSGARHLESEKLA QILAARELEIKQIPSDGHCMYGALEDQLREQDCALTVASLRRQTAEYMQTHSDDFLPFLTNPSTGDMYTP EEFGKYCDDIVNTAAWGGQLELRALSHILQTPIEILQADAPPIIVGEEYPRNPLVLVYMRHAYGLGEHYN SVTRLVNSATENC SGPTRRRLLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC031474
ORF Size:	882 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC031474.1
RefSeq Size:	2598 bp
RefSeq ORF:	884 bp
Locus ID:	72201
Cytogenetics:	4 A1
MW:	33.8 kDa