

Product datasheet for **MR204053L3V**

Rel2 (NM_153793) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rel2 (NM_153793) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rel2
Synonyms:	4631403P03Rik; ependolin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153793
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204053).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153793.1 , NP_722488.1
RefSeq Size:	1783 bp
RefSeq ORF:	912 bp
Locus ID:	225392
UniProt ID:	Q8BRJ3
Cytogenetics:	18 B3
Gene Summary:	Induces activation of MAPK14/p38 cascade, when overexpressed. Induces apoptosis, when overexpressed.[UniProtKB/Swiss-Prot Function]



[View online »](#)