

Product datasheet for **MR204055L3V**

2310057J18Rik (BC021783) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: 2310057J18Rik

Synonyms: 2210414P04Rik; Leg1

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: BC021783

ORF Size: 882 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR204055).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC021783](#), [AAH21783](#)

RefSeq Size: 1250 bp

RefSeq ORF: 884 bp

Locus ID: 67719

Cytogenetics: 10 A4

Gene Summary: May be involved in early liver development.[UniProtKB/Swiss-Prot Function]

