

## Product datasheet for **MR204062L4V**

### Utp15 (BC034372) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Utp15 (BC034372) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Utp15
Synonyms:	AW544865
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC034372
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204062).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC034372</a> , <a href="#">AAH34372</a>
RefSeq Size:	2209 bp
RefSeq ORF:	887 bp
Locus ID:	105372
Cytogenetics:	13 D1
Gene Summary:	Ribosome biogenesis factor. Involved in nucleolar processing of pre-18S ribosomal RNA. Required for optimal pre-ribosomal RNA transcription by RNA polymerase I.[UniProtKB/Swiss-Prot Function]


[View online »](#)