

Product datasheet for **MR204068A1V**

Mouse Asb17 (NM_025758) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Asb17 (NM_025758) AAV Particle
Tag:	Myc-DDK
Symbol:	Asb17
Synonyms:	4933436O18Rik
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204068 representing NM_025758 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATAACTCTTCTAAATTATGCCGGAAGACGTCTTTTCCAAGAAGCAATATATTCTGTAACTTGTTG
ACAAGATAGTTAAGCGGCCCTCTCTGCAGTTTCTGGCCAATGGGGATACCACTGCTATGAACCCAGGAT
TTACAGAACCCTGGCAAAATCCTGAGGTATGTTGACTTGGACGGCTTTGACATACTCCTCACGGACTAT
ATTGCTTTTGTGAAAAGTCAGGACACCGTTTGAACCTCACTTTAACCTTGAGTTTACTGAAATATGCG
TGAATACCATCTGTACTGGGTTTTGCCAGGAAAGGTAATCCTGACTTCGTGGAAGTCTTCTCAAGAA
GACGAAGGACTATGTCCAAGACAGAAGCTGCAGCCTGGCGCTGATATGGAGAACCCTTCACACCTGTGTAC
TGCCCCAGCCCCCTGAGTGGCATCACACCTCTACTCTACGTGGCTCAGACAAGACAGTCAAATATCTTAA
AAATTCTCCTGCAGTATGGAATCCTAGAAAGAGAAAAAACCTATCAACATTGTTCTGACAATACTACT
TTACCCTTCGAGAGTGAGAATAATGGTTGACCAGAGTTGATTGACATTCAAGAAGATGCCAAGACATGT
TTAATGCTATGTTCCAGAGTGCTTTCTACGATCTCAGTCAGGGAGATAGAGACACAGCTGAGCTTAGGAC
GACGCCCAATTATTCAAATTGGTTGGACTACATCCCGCCAACAAGATAACAAGGATCCATGTGAACCTCGT
CCACCTTTCAGATAAACCATCAGGACCAACTGCTGGCCAACAATATGCTCCCAATGGAATATTTTCC
CTTCTAATTCCTACTCGTTTACAAAACCTCCTGAATTTAGAAAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204068 representing NM_025758 Red=Cloning site Green=Tags(s) MNSSKLCRKTSFPRSNI FCNLVDKIVKRPSLQFLGQWGYHCYEPRIYRTLAKILRYVDLDGFDILLTDY IAFVEKSGHRFELNFNLEFT EICVNTILYWVFARKGNPDFVELLLKTKDYVQDRSCSLALIWRFTPVY CPSPLSGITPLLYVAQTRQSNILKILLQYGILEREKNPINIVLTILLYPSRVIMVDHELIDIQEDAKTC LMLCSRVLSTISVREIETQLSLGRRPIIQNWLDYIPTRYKDPCELVHLCRITIRTQLLANMPLPNGIFS LLIPTRLQNFLNLES TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_025758
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_025758.4, NP_080034.2</u>
RefSeq Size:	1048 bp
RefSeq ORF:	888 bp
Locus ID:	66772
UniProt ID:	<u>Q8VHP9</u>
Cytogenetics:	3 H3
MW:	34.5 kDa