

Product datasheet for **MR204070A1V**

Mouse Sfrp2 (NM_009144) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Sfrp2 (NM_009144) AAV Particle
Tag: Myc-DDK
Symbol: Sfrp2
Synonyms: AI851596; Sdf5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR204070 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGCGGGGCCCTGCCTCGCTGCTGCTAGTCCTCGCCTCGCACTGCTGCCTGGGCTCGGCGCGTG
GGCTCTTCCTCTTCGGCCAGCCCGACTTCTCTACAAGCGCAGCAACTGCAAGCCCATCCCCGCCAACCT
GCAGCTGTGCCACGGCATCGAGTACCAGAACATGCGGCTGCCAACCTGCTGGGCCACGAGACCATGAAG
GAGGTGCTGGAGCAGGCGGGCGCCTGGATTCCGCTGGTCATGAAGCAGTGCCACCCGGACACCAAGAAGT
TCCTGTGCTCGCTCTTCGCCCTGTCTGTCTCGACGACCTAGATGAGACCATCCAGCCGTGTCACCTCGCT
CTGCGTGCAGGTGAAGGACCGCTGCGCCCCGGTCATGTCGCTTCGGCTTCCCTGGCCAGACATGCTG
GAGTGCGACCGTTTCCCGCAGGACAACGACCTCTGCATCCCCCTCGCTAGTAGCGACCACTCCTGCCGG
CCACAGAGGAAGCTCCCAAGGTGTGAAGCCTGCAAAACCAAGAATGAGGACGACAACGACATCATGGA
AACCTTTGTAATAATGACTTCGCACTGAAAATCAAAGTGAAGGAGATAACGTACATCAACAGAGACACC
AAGATCATCCTGGAGACAAAGAGCAAGACCATTTACAAGCTGAACGGCGTGTCCGAAAGGGACCTGAAGA
AATCCGTGCTGTGGCTCAAAGACAGCCTGCAGTGCACCTGTGAGGAGATGAACGACATCAACGCTCCGTA
TCTGGTCATGGGACAGAAGCAGGCGGGCAGCTGGTGATCACCTCCGTGAAACGGTGGCAGAAGGGCCAG
AGAGAGTTCAAGCGCATCTCCCGCAGCATCCGCAAGCTGCAATGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204070 protein sequence Red=Cloning site Green=Tags(s) MPRGPASLLLLVLASHCCLGSARGLFLFGQPDFSYKRSNCKPIPANLQLCHGIEYQNMRLPNLLGHETMK EVLEQAGAWIPLVMKQCHPDTKKFLCSLFAPVCLDDLDETIQCHSLCVQVKDRCAPVMSAFGFPWPDML ECDRFPQDNDLCIPLASSDHLLPATEEAPKVCEACKTKNEDDNDIMETLCKNDFALKIKVKEITYINRDT KIILETKSKTIYKLVNGVSRDLKKSVLWLDLSLQCTCEEMNDINAPYLVMGQKQGELVITSVKRWQKGQ REFKRISRSIRKLQC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_009144
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_009144.1</u>
RefSeq Size:	2001 bp
RefSeq ORF:	888 bp
Locus ID:	20319
UniProt ID:	<u>P97299</u>
Cytogenetics:	3 37.37 cM
MW:	33.5 kDa