

Product datasheet for **MR204082A1V**

Mouse Blvra (NM_026678) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Blvra (NM_026678) AAV Particle
Tag:	Myc-DDK
Symbol:	Blvra
Synonyms:	0610006A11Rik; 2500001N03Rik; Blvr
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204082 representing NM_026678 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTACTGAGCCAAAGAGGAAATTTGGTGTGGTAGTGGTTGGTGGTGGCAGAGCTGGCTCTGTGAGGA
TAAGGGACTTGAAGGATCCCACTCTTCAGCATTCTTAAACCTGATTGGATATGTGTCCAGACGAGAACT
TGGGAGCCTTGATAATGTACGGCAGATTTCTTTGGAAGATGCTCTCCGAGCCAGAGGTTGATGTGCC
TATATTTGCACTGAGAGTTCAGCCACGAAGACTATATAAGGCAGTTTCTTCAGGCTGGCAAGCATGTCC
TCGTGGAATACCCCATGGCATTGTCAATTTGCGGCAGCGCAGGAGCTGTGGGAGCTGGCTGCACAGAAAGG
GAGAGTCCTGCATGAGGAGCACATCGAACTCTTGATGGAAGAGTTCTGAATTCCTGAAAAGAGAAGTGGCC
GGGAAGGAGCTGCTAAAAGGATCTCTTCGCTTCACAGCCAGCCCTGGAAGAAGAGAAATTTGGCTTCC
CTGCGTTCAGTGGCATTCTCGCCTGACCTGGCTGGTCTCCCTCTTTGGGAGCTTTCTCTATTCTGC
CACCATGGAACCGAAAAGAAGATCAATATATGAAAATGACCGTGCAGCTCGAGACCCAGAACAAAAGT
CCATTATCATGGATTGAAGAGAAAGGGCTGGCTTAAAAGGAACAGACATATAAGCATCCACTTCAAAT
CTGGGTCCTTGAGGAAGTGCCAAATGTAGGAGTCAATAAGAACATTTCTGAAAGATCAGGATATATT
TATCCAGAAGCTCTTGGGCCAGGTCTCTGCAGAGGACCTGGCTGCTGAGAAGAAGCGCATCTGCATTGC
CTGGAGCTGGCCAGCGACATCCAGAGGCTGTGCCACCGGAAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204082 representing BC052146 Red=Cloning site Green=Tags(s) MSTEPKRKFGVVVVGVRAGSVRIRDSKDPHSSAFLNLIGYVSRRELGSLDNVRQISLEDALRSQEVDVA YICTESSSHEDYIRQLQAGKHVLVEYPMALSFAAAQELWELAAQKGRVLHEEHIELLMEEFEFLKREVA GKELLKGSRLFTASPLEEEKFGFPAFSGISRLTWLVSLFGELSLISATMENRKEDQYMKMTVQLETQNK PLSWIEEKGPGLKRNHRHSIHFKSGSLEEVPNVGVNKNIFLKDQDIFIQKLLGQVSAEDLAAEKKRILHC LELASDIQRLCHRKQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026678
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026678.4</u>
RefSeq Size:	1394 bp
RefSeq ORF:	888 bp
Locus ID:	109778
UniProt ID:	<u>Q9CY64</u>
Cytogenetics:	2 61.76 cM
MW:	33.5 kDa