

Product datasheet for **MR204099L4V**

Arrdc4 (NM_025549) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arrdc4 (NM_025549) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arrdc4
Synonyms:	2410003C09Rik; AV216361
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025549
ORF Size:	891 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204099).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025549.3 , NP_079825.2
RefSeq Size:	1097 bp
RefSeq ORF:	891 bp
Locus ID:	66412
UniProt ID:	A0A0B4J1F4
Cytogenetics:	7 C



[View online »](#)

Gene Summary:

Functions as an adapter recruiting ubiquitin-protein ligases to their specific substrates (PubMed:27462458). Plays a role in endocytosis of activated G protein-coupled receptors (GPCRs) (By similarity). Through an ubiquitination-dependent mechanism plays also a role in the incorporation of SLC11A2 into extracellular vesicles (PubMed:27462458). May play a role in glucose uptake (By similarity).[UniProtKB/Swiss-Prot Function]