

## Product datasheet for **MR204120L4V**

### Abhd10 (NM\_172511) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Abhd10 (NM_172511) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abhd10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172511
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204120).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172511.2</a>
RefSeq Size:	2860 bp
RefSeq ORF:	894 bp
Locus ID:	213012
UniProt ID:	<a href="#">Q6PE15</a>
Cytogenetics:	16 B5
Gene Summary:	This gene encodes a mitochondrially-localized enzyme that acts as a hydrolase. The encoded protein removes glucuronide from mycophenolic acid acyl-glucuronide. There are two pseudogenes for this gene on chromosome 8. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013]



[View online »](#)