

Product datasheet for **MR204132**

AI118078 (BC068128) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	AI118078
Synonyms:	9630029F15; AI118078; HVRP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR204132 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGCTATTTAGCCACACCGTGCTTTTCAAGTCTGACTGGCTGAAGCCGTGCTGTGGGAAGAGAGCAG
CCGTATGGCAGGTATTTTGTCTAGTGCAAGTCTCAACAGTTTCTGGTAGCCTGTGTAATATTGGTGGT
GATCCTCCTGACTCTGGAGCTTCTCATAGATACAAAGCTTCTCCAGTTTCCAATGCTTCCAGTTTGCT
GGTGTCATTCACTGGATCAGTCTGGTCATTCTCTGTGTTCTTCTCAGAGACTGCTCCTACGGATCGTGG
TACTGGGGATCTGGGATTACATCGAAAACAAAATAGAGGTGTTTGACGGGGCTGTGATCATCTGTCTT
GGCCCCGATGGTGGCGTCCACTGTGGCTAACGGACCCAGGAGCCCTGGGATGCCATCAGCCTCATCATC
ATGTTCCGAATCTGGCGGTGAAGAGGGTCATTGATGCCTATGCTCTGCCAGTCAAGTTGGAGATGGAGA
TGGTCAACCAGCAGTATGAGAAGGCCAAGGCCATCCAAGATGAGCAGCTGGAAGACTGACGCAATCTG
TCAGGAGCAAGGGTTTGAATCCGGCAGCTGCGTGCGCACCTGGCACAGCAGGACCTGGATCTGGCAGCC
GAGCGGGAGGCGCGCTGCAGGCCCCACACGTGCTCAGCCAGCCACGCAGCCGCTACAAGGTCGTAGAGG
CTGGCACGTGGGCCGAGGAGACAGCAGCCGAGAGCATCGTGAAGAGCTGAGGCCCTCTCAAGAAGCCAC
AGTGAAGATGACATGAACAGCTACATCAGCCAATACTACAATGGGCCAGCAGTGGTAAGTCTGGGTTG
GGACTGTTTGCATGCACTATGCAGCCCCAGGACACTTGTGACCCAGCAAAACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR204132 protein sequence
 Red=Cloning site Green=Tags(s)


MMLFQPHRAFRSDWLKPCCGKRAAVWQVFLLSASLNSFLVACVILVVILLTLELLIDTKLLQFSNAFQFA
 GVIHWISLVLSVFFSETVLRIVVLGIWDYIENKIEVFDGAVIILSLAPMASTVANGPRSPWDAISLII
 MFRIWRVKRVIDAYVLPVKLEMEMVTQQYEKAKAIQDEQLERLTQICQEQGFIRQLRAHLAQDDLAA
 EREAALQAPHVLSQPRSRYSKVEAGTWAEEETAESIVEELRPSQEATVKDDMNSYISQYYNGPSSGKSGL
 GLFACTMQPQDTCPSKP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

BC068128

ORF Size:

894 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC068128](#), [AAH68128](#)

RefSeq Size: 4432 bp

RefSeq ORF: 896 bp


Locus ID: 244886

Cytogenetics: 9 B

MW: 33.6 kDa

Gene Summary: Voltage-sensor protein present on the post-synaptic side of glutamatergic mossy fibers and granule cells in the cerebellum. Despite the presence of a voltage-sensor segment, does not form a functional ion channel and its precise role remains unclear. Undergoes both rapid and slow structural rearrangements in response to changes in voltage. Contains a zinc-binding site that can regulate the slow conformational transition.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204132