

## Product datasheet for **MR204133L4V**

### Hsd17b11 (NM\_053262) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsd17b11 (NM_053262) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hsd17b11
Synonyms:	Dhrs8; Pan1b; retSDR2; SDR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053262
ORF Size:	897 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204133).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053262.3</a>
RefSeq Size:	1699 bp
RefSeq ORF:	897 bp
Locus ID:	114664
UniProt ID:	<a href="#">Q9EQ06</a>
Cytogenetics:	5 E5



[View online »](#)

**Gene Summary:**

Can convert androstan-3-alpha,17-beta-diol (3-alpha-diol) to androsterone in vitro, suggesting that it may participate in androgen metabolism during steroidogenesis. May act by metabolizing compounds that stimulate steroid synthesis and/or by generating metabolites that inhibit it. Has no activity toward DHEA (dehydroepiandrosterone), or A-dione (4-androste-3,17-dione), and only a slight activity toward testosterone to A-dione. [UniProtKB/Swiss-Prot Function]