

Product datasheet for **MR204136L4V**

Fam175a (Abraxas1) (NM_172405) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam175a (Abraxas1) (NM_172405) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abraxas1
Synonyms:	3830405G04Rik; 5630400M01Rik; AI506069; AL024423; AV118690; Ccdc98
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172405
ORF Size:	897 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204136).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172405.3 , NP_765993.1
RefSeq Size:	2499 bp
RefSeq ORF:	1224 bp
Locus ID:	70681
UniProt ID:	Q8BPZ8
Cytogenetics:	5 E4



[View online »](#)

Gene Summary:

Involved in DNA damage response and double-strand break (DSB) repair. Component of the BRCA1-A complex, acting as a central scaffold protein that assembles the various components of the complex and mediates the recruitment of BRCA1. The BRCA1-A complex specifically recognizes 'Lys-63'-linked ubiquitinated histones H2A and H2AX at DNA lesion sites, leading to target the BRCA1-BARD1 heterodimer to sites of DNA damage at DSBs. This complex also possesses deubiquitinase activity that specifically removes 'Lys-63'-linked ubiquitin on histones H2A and H2AX.[UniProtKB/Swiss-Prot Function]