

Product datasheet for **MR204149**

Khk (NM_008439) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Khk (NM_008439) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Khk |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >MR204149 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGAGAAGCAGATCCTGTGCGTGGGGCTGGTGGTCTGGACATCATCAATGTGGTGGACAAATACC
CAGAGGAAGACACGGATCGCAGGTGCCTGTCCCAGAGATGGCAGCGTGGAGGCAACGCATCCAACCTCTG
CACTGTCCTTTCCTTGCTTGGAGCCCGCTGTGCCTTCAATGGCTCTTTGGCCCTGGCCACGTTGCCGAC
TTCCTGGTGGCTGACTTCAGGCAGAGGGCGTGGATGTGTCTCAAGTGACTTGGCAGAGCCAGGGAGATA
CCCCTTGTCTTGTGTCATCGTCAACAACCTCAATGGCTCCCGTACCATTATACTCTACGACACGAACT
GCCAGATGTGTCTGCTAAGGACTTTGAGAAGGTCGATCTGACCCGGTTCAAGTGGATCCACATTGAGGGC
CGGAATGCATCGGAACAGGTGAAGATGCTGCAGCGGATAGAGGAGCACAATGCCAAGCAGCCTCTGCCAC
AGAAGGTCGGGTGTCGGTGGAGATAGAGAAGCCCGTGAGGAGCTCTTCCAGTTGTTTAGCTATGGTGA
GGTGGTGTGTCAGCAAAGATGTGGCCAAGCACCTGGGGTTCAGTCAGCAGTGGAGGCCCTGAGGGGC
TTGTACAGTCGAGTGAAGAAAGGGGCTACGCTTGTCTGTGCTGGGCTGAGGAGGGTGCCGATGCCCTGG
GCCCGATGGTCAGCTGCTCCACTCAGATGCCTTCCCACCGCCCCGAGTAGTAGACTCTTGGGGCTGG
AGACACCTTCAATGCCTCTGTCTTTCAGCCTCTCGAAGGGAACAGCATGCAAGAGGCCCTGAGATTC
GGTGCCAGGTGGCTGGCAAGAAGTGTGGCTTGCAGGGGTTTATGGCATTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204149 protein sequence
Red=Cloning site Green=Tags(s)

MEEKQILCVGLVVLDIINVVDKYPEEDDRRCLSQRWQRGGNASNSCTVLSLLGARCAFMGSLAPGHVAD
 FLVADFRQRGVDVSQVTWQSQGDTPCSCCIVNNSNGSRTIILYDTNLPDVSADFEKVDLTRFKWIHIEG
 RNASEQVKMLQRIEEHNAKQPLPQKVRVSVIEIEKPREELFQLFSYGEVVFVSKDVAKHLGFQSAVEALRG
 LYSRVKKGATLVCAWAEEGADALGPDGQLLHSDAFPPPRVVDLTLGAGDTFNASVIFSLSKGNMQEALRF
 GCQVAGKKCGLQGFDFGIV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_008439

ORF Size: 897 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008439.4](#)

RefSeq Size: 1117 bp

RefSeq ORF: 897 bp

Locus ID: 16548


UniProt ID: [P97328](#)

Cytogenetics: 5 16.9 cM

MW: 32.7 kDa

Gene Summary: Catalyzes the phosphorylation of the ketose sugar fructose to fructose-1-phosphate. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204149