

Product datasheet for **MR204150L4V**

Slc25a4 (BC026925) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a4 (BC026925) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a4
Synonyms:	Ant1; AU019225
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC026925
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204150).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC026925 , AAH26925
RefSeq Size:	1070 bp
RefSeq ORF:	896 bp
Locus ID:	11739
Cytogenetics:	8 26.22 cM
Gene Summary:	Involved in mitochondrial ADP/ATP transport. Catalyzes the exchange of cytoplasmic ADP with mitochondrial ATP across the mitochondrial inner membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)